

# Last round of DCAM specimens set for launch


## Space Shuttle-Mir Science Program (STS-89/91)

Experiments with the Diffusion-controlled Crystallization Apparatus for Microgravity (DCAM) aboard the Mir space station will conclude with a payload of six trays of cells to be carried by STS-89 in January 1998.

To date, with stays lasting as long as 6-months aboard Mir, DCAM has yielded dramatic results. Highlights include numerous large, spectacular crystals of the nucleosome core particle (shown at right), which regulates genetic activities in the nucleus of a cell. Another striking result was the growth of the largest T7 RNA crystal ever produced (0.7 x 0.8 x 2.0 mm in size). DCAM cells carrying triglycine sulfate (TGS) also yielded large crystals. TGS has been grown in space by other experiments, so these results will help in gauging the comparative effectiveness of various microgravity processes.

DCAM's first flight, the second U.S. Microgravity Laboratory (USML-2) in 1995, verified the scientific validity of this approach to growing protein crystals. A full complement of six trays comprising 162 DCAM cells was carried to Mir in March 1996 and was returned to Earth on STS-79 which also installed a second set of DCAM trays. The STS-81 mission in January 1997 retrieved the second set of trays and installed a third set, retrieved by STS-85 in May 1997, as scientists continued to refine the mixtures and details used in this promising method. STS-91 will retrieve this DCAM set in May 1998.

Proteins are important, complex biological molecules which serve a variety of functions in all living organisms. Determining their three-dimensional structure will lead to a greater understanding of how they function in living organisms. Many proteins can be crystallized and their molecular structures determined through analysis of those crystals by X-ray diffraction. Unfortunately, some



Crystal of nucleosome core particle.


crystals grown in the 1-g environment of Earth have internal defects that limit or impair such analyses. As demonstrated on Space Shuttle missions since 1985, some protein crystals grown in space are larger, and more highly ordered, than the Earth-grown counterparts.

DCAM was developed at Marshall Space Flight Center to grow protein crystals by a special diffusion process. The principal investigator is Dr. Daniel Carter of New Century Pharmaceuticals.

### ► Diffusion-controlled Crystallization Apparatus for Microgravity (DCAM)

DCAM grows crystals by the dialysis and liquid / liquid diffusion methods. In both methods, protein crystal growth is induced by changing conditions in the solution containing the protein. In dialysis, a semipermeable membrane retains the protein solution in one compartment, while allowing molecules of precipitant to pass freely through the membrane from an adjacent compartment. As the precipitant concentration increases within the protein compartment, crystallization begins.

In liquid-liquid diffusion, a protein solution and a precipitant solution are layered in a container and allowed to diffuse into each other. This leads to conditions which may induce crystallization of the protein. Liquid-liquid diffusion is difficult on Earth




Biotechnology Program Office



Protein Crystal Growth Program


#### Hardware

**DCAM** (Diffusion-controlled Crystallization Apparatus for Microgravity)



Marshall Space Flight Center  
Huntsville, Alabama

The DCAM has no mechanical system. Diffusion starts on Earth as soon as the chambers are filled. However, the rate is so slow that no appreciable change occurs before the samples reach orbit one or two or even several days later. This also allows protein samples to stay aboard the shuttle in case of a launch delay. In other hardware, many samples must be replaced in the event of a postponement. Such an apparatus is ideally suited for long



STS-89 will carry 162 DCAM units mounted in 3 x 9 arrays on six trays stored in a locker in the Shuttle middeck (the same as the array now aboard Mir). Upon arrival at Mir, the DCAMs will be transferred to one of Mir's modules and mounted in a quiet area where crystallization will take place. Three of the six trays will be mounted so they can be photographed as crystals form. After the retrieval and return to Earth by STS-92 in May 1998, the samples will be analyzed.

## Outer surface glycoprotein of the hyperthermophile *Methanothermus fervidus*

**Ferrochelata** is important to biomedical and biochemical applications. *Investigators:* Dr. B.C. Wang and Dr. Harry Dailey University of Georgia.

[illegible]